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PROBLEM	CAUSE	SOLUTION
Loud noise	<ol style="list-style-type: none"> 1. Over speed rotor 2. Device is tilting 3. Unstable surface 4. Loose mechanical components 5. Spindle motor or motor are bent 	<ol style="list-style-type: none"> 1. Be sure not to exceed the maximum rated speed of the rotor (Table 1) 2. Put the device in a horizontal and stable place 3. Put the device in a horizontal and stable place 4. Screw the rotor 5. Contact with Nahita technical assistance 6. Check the spindle and the rotor after screw 7. Contact with Nahita technical assistance
Unstable turn speed	Fault of control circuit	Contact with Nahita technical
Motor failure (smell of burnt, etc.)	<ol style="list-style-type: none"> 1. Fault of control circuit 2. Electric control or motor are burnt 	Contact with Nahita technical

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, does not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

SYMBOL	PROBLEM	SOLUTION
E1	Imbalance. Centrifuge stop run because of over-vibration	<ol style="list-style-type: none"> 1. Re-weight error allowed $\leq 1.5g$ 2. The device is placed unbalanced and the stress is unequal, please adjust the device to make the stress equal 3. Clean any dirty on the surface of the rotor or inside of the tube holes 4. If the spindle motor is bent, please contact with Nahita technical assistance
E3	Protección de tapa	<ol style="list-style-type: none"> 1. Vuelva a cerrar correctamente la tapa 2. Si la tapa está rota contacte el Servicio Técnico Nahita
E7	Fallo de motor	Contacte el Servicio Técnico Nahita
E8	Selección errónea de parámetros	Reseteo los parámetros de acuerdo a la Tabla 1

PROBLEM	CAUSE	SOLUTION
No display	<ol style="list-style-type: none"> 1. Wrong connection 2. Electrical interference 3. The fuse is blown 4. Bad connection or motherboard failure 	<ol style="list-style-type: none"> 1. Check whether the power outlet and the connection is ok. 2. Switch off and switch on after 1 min 3. Change the fuse 4. Contact with Nahita technical assistance
Display is ok but device doesn't run	<ol style="list-style-type: none"> 1. Door is not good close 2. Door is not open after centrifuge 3. The motor doesn't work after pressing Start key 	<ol style="list-style-type: none"> 1. Close correctly the door 2. Open the door, re-close it correctly and try again 3. Contact with Nahita technical assistance
Door can't open	Electrical fault	<ol style="list-style-type: none"> 1. Open de door manually 2. Contact with Nahita technical assistance
Unbalance alarm	<ol style="list-style-type: none"> 1. Samples are not placed in a balanced position 2. Tube adapter and tube holder can't swing freely 3. Dirty on the surface of the rotor or inside of the tube holes 	<ol style="list-style-type: none"> 1. Place again the tubes in a symmetric position 2. Check adapters position 3. Clean dirty pieces

Thank you for choosing this equipment. We sincerely wish that you enjoy your Centrifuge Nahita Multibas 2750. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Centrifuge Nahita Multibas 2750, so it should be available for all the users of this equipment.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ You should carefully handle the centrifuge avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the centrifuge to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the centrifuge.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es.
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the centrifuge warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the centrifuge warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

SPANISH	2-15
ENGLISH.....	16-28
FRENCH.....	29-42

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT.....	17
2. DESCRIPTION	18
3. TECHNICAL SPECIFICATIONS	19
4. THEORETICAL BASIS	19
5. INSTALLATION / SETTING UP.....	19
6. MAINTENANCE AND CLEANING.....	24
7. ACCESSORIES	25
8. TROUBLESHOOTING.....	26


1. USES OF THE INSTRUMENT

Nahita multibas 2750 centrifuge is a tabletop, silence and vibration-free construction. The centrifuge is designed for a low-speed use and provided with a swing rotor and different tube adaptors for it. It offers the option to adapt a second special rotor for microtiter plates.


9 different operating procedures for choose with continuous control of operating setting through its wide LCD display. The centrifuge includes electromechanical lock in the door cover, alarm, and automatic stop for over speed or rotor imbalance and temperature control of the microprocessor as a security system.

This centrifuge can be widely used in fields of clinical medical, biology, chemistry, genetic engineering, immunology, etc. due to its simple design and the quality of its components. It is also able to be operated for long periods of time or to maintain a high rate of cycles per day.

Cabecal oscilante

Cabecal placas microtiter

8. TROUBLESHOOTING

The table below shows the alarm information of the machine (2.13) and the reason as well as the way to solve; you can rule it out according to the table. If users still can't solve the problem or the alarm information is not included, please contact maintainer from Nahita Technical Assistance immediatel

SYMBOL	PROBLEM	SOLUTION
E2	Over speed. It stops to work because of detection of over speed of rotor	1. Problem of microcomputer control system, please contact with Nahita technical assistance 2.Problem of speed sensor, please contact with Nahita technical assistance

Cleaning

Working at high speed exist the possibility that the test-tubes become broken. If this happens, stainless steel cavity and adapters must be cleaned and disinfected. The broken pieces should also be removed. If you do not perform this operation will be damage the different components of the centrifuge.

Please turn off the power switch and unplug the power cord before the doing cleaning or maintenance of the centrifuge.

Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminum, coatings, etc. as they damage the centrifuge and produce an early ageing of the equipment.

Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

7. ACCESSORIES


OPTIONAL ADAPTERS NOT INCLUDED.

Code	Description
92750001	Adapter 4x100 mL P/Multibas
92750002	Adapter 4x50 mL P/Multibas
92750003	Adapter 8x50 mL P/Multibas
92750004	Adapter 24x15 mL P/Multibas
92750005	Adapter 32x15 mL P/Multibas
92750007	Adapter 48x7 mL P/Multibas
92750008	Adapter 72x7 mL P/Multibas
92750009	4 Microtiter plates rotor P/Multibas
92750010	Adapter 12x15 mL P/Multibas
92750011	Adapter 4x250 mL P/Multibas

Rotor code	Capacity	RPM,	Máx. FCR	Tubes (ØxL)
0001 (swing)	4 x 50 mL	5000	4980	Falcon (30x117 mm)
0002 (swing)	4 x 100 mL	5000	4600	42x100 mm
0003 (swing)	8 x 50 mL	4000	3040	Falcon (30x117 mm)
	24 x 15 mL	4000	3040	Falcon (16.5x120 mm)
	32 x 15mL	4000	3040	Falcon (16.5x120 mm)
	48 x 7 mL	4000	3040	Blood extraction tube (13x100 mm)
0004 (swing)	72 x 7 mL	4000	3000	Blood extraction tube (13x100 mm)
0005 (swing)	2x2x96 Pocillos	3500	2120	Microtiter plate
0006 (swing)	12 x 15 mL	5000	3500	Centrifuge essay tube
0007 (swing)	4 x 250 mL	4000	3010	Centrifuge essay tube


2. DESCRIPTION

- 1.1 Door cover
- 1.2 Fuselage sell
- 1.3 Suction cup legs
- 1.4 Switch
- 1.5 Cord for emergency lock
- 1.6 Front panel
- 1.7 Front air intake
- 1.8 Control panel (Display /keypad)
- 1.9 Observation hole
- 1.10 Power outlet

Front panel

- 2.1 Open door key
- 2.2 Stop key
- 2.3 Start/moment centrifugal key
- 2.4 Setting knob
- 2.5 Acceleration/deceleration settings
- 2.6 Program number
- 2.7 Temperature display of driver board (security control)
- 2.8 Set the centrifugal time
- 2.9 Set speed (RPM) and relative centrifugal force (RCF)
- 2.10 Showing real-time RPM and RCF while running
- 2.11 Software version
- 2.12 Showing centrifugal time count-down
- 2.13 Sign of fault/display error
- 2.14 Rotor number
- 2.15 Status display

3. TECHNICAL SPECIFICATIONS

Code	52750000
Rotor	Swing
Speed	100-5000 rpm (100 to 100)
Max. RCF	4980 g (Table 1)
Timer	1'-1h 59' 59"
Noise	<65dB
Dimensions (L x W x H)	60 x 47 x 35 cm
Weight	30 kg
Power supply	AC230-240 V, 50/60 Hz
Fuse	10 A


4. THEORETICAL BASIS

Relative centrifugal force is based in the particle movement. RCF is usually thousands of times of gravity (g). It's unit items to measure the efficiency of variety of instrument separation or precipitation. The calculation of centrifugal force has something to do with centrifugal speed and centrifugal radius. It's based on the following equation:

$$F.C.R.= 1.118 \cdot 10^{-5} \cdot r \cdot n^2$$

r The centrifugal radius (cm).

n Centrifugal speed (rpm).

5. INSTALLATION / SETTING UP

Preliminary inspection

- ◆ Unwrap the balance, take off the involving plastic and take off the poliespan protection in which it comes fitted.
- ◆ Please make sure that it does not present any damage because of the shipment. If it does please communicate to your transport agent or suministrator immediately, so that they can make the claims in the correct time limit.

We will only accept any equipment return within 15 days after delivery and provided it comes in its original package and with all accessories and documents supplied

Please check that all the accessories are enclosed with the equipment:

- Rotor
- ◆ - 6 mm hexagon wrench
- Power cord

- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the centrifuge by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the centrifuge (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

6. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ The centrifuge tubes can load the samples at different time but it must be symmetric loading of samples (the weight of error allowed $\leq 1,5$ g). Start the machine with asymmetric loading is not allowed.
- ◆ Usually, rotor should be clean once a week. For separating salts or corrosive samples, please clean immediately after using. If the sample spills out or drip in the rotor, please blot it and do partial cleaning.
- ◆ Please ensure that proper rotor are installed tighten before the running of the centrifuge.
- ◆ Don't open the door manually if the centrifuge is still running (rotor is turning) or the rotors don't stop to turn
- ◆ Choose proper container resistant to PCF used for placing the samples.
- ◆ Check whether the centrifugal container and parts are damaged or not when it's under maximum allowed load and maximum allowed speed. If there is anything damaged, please replace in time.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The centrifuge is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the centrifuge by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the centrifuge (burns, hurts...) or damages in nearby equipment.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Note: If during the process of setting parameters, if there is alarm caused by machine fault or incorrect setting on the parameters, press the stop key (2.2) to cancel the alarm.

◆ Press the start key (2.3). The machine will begin to run and goes to the setting speed. You can see the countdown time (2.12) and the speed and RCF (2.10) while working.

When the centrifuge starts running the electromechanical lock is blocking the cover door, never force its opening in this situation.

When the rotating speed of the rotor is over 250 rpm of setting speed, the machine will give an alarm; when the rotating speed is over 300 rpm of maximum rated speed of the rotor, the rotor will stop to rotating automatically and you can open the door until the rotor stops totally.

Note: Press stop key to stop the alarm (2.2)

◆ If you need stop the program at any time, press stop key (2.2).

Time setting will reduce from set value to zero, when it shows zero, the centrifuge will stop working. And at this time, the speed will reduce from set value to zero. When the speed is zero, the door cover will open automatically and the machine will buzz, please press stop key (2.3) to stop it and centrifugal process is completed. It is recommended switch off the device from the main knob (1.4) to avoid any accident.

If you need short time centrifugal, please press start key (2.3) for long time, and the speed will increase gradually. It will stop when you release and the maximum speed is the setting speed of the rotor.

Once the device has stabilized the speed of work you can modify the parameters setting (the changes will set after 2-3 seconds)

Please note the using life of the centrifugal container. Check whether the centrifugal container (plastic, glass) is damaged or not when it's under maximum allowed load and maximum allowed speed. If there's anything damaged, please replace in time

Security

◆ The centrifuge must be used by previously qualified staff that knows how the equipment works thanks to the user manual.

◆ Is heavy equipment. Please lift the machine from the left and right sides balances. Please move it vertically and don't swag it.

◆ You should put the centrifuge in a horizontal plane stable table, having a safety area of at least 30 cm per side.

◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).

◆ Do not place the centrifuge near any warm supply (burners, blowlamps, etc.), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.

◆ Do not use following materials in the centrifuge: Flammable and explosive materials, strong chemical materials and toxic or radioactive.

◆ During its functioning dangerous materials such as flammable or pathological substances must be out of the safety area.

◆ When you are not using the centrifuge for a long period of time please make sure it is unplugged in order to avoid possible accidents.

- User guide

Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

Is heavy equipment, please lift the machine from the left and right sides balances. Please move it vertically and don't swag it.

◆ Please put the centrifuge on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the centrifuge near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.

◆ Avoid inflammable or toxic substances in the working area.

◆ The centrifuge is supplied with a Schuko standard wire

◆ Please insert the wire that feeds the AC electric current in the base of current 230-240 V, 50-60Hz $\pm 10\%$ provided with earth wire and to the other end to the centrifuge connector.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 230-240 V 50-60Hz $\pm 10\%$.

◆ If you are not using the device for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

◆ The environmental temperature for the room keeping the centrifuge should be between 5 °C – 40 °C, and the environmental humidity should be <80%. Please keep the environmental clean.

In order to prevent damage to the spindle, please remove the rotor out before moving the centrifuge.

Don't cut the electricity when the rotor is running, or it will damage the control circuit.

Switch on

◆ With the door closed, switch on (1.4) and wait until self-examination (red light on control panel). Display shows main work panel.

◆ Press (2.1) and the door cover will eject automatically to a certain height. Check that the inner space is clean and remove any external object.

Please prevent your head or anything else from the upward side of the door cover to avoid dangers.

If there is power failure suddenly or there are troubles when the rotor is rotating, and you can't open the door by pressing keys, please open the door manually (1.5). Proceed as follow:

- Check through the observation hole (1.9) that the rotor is stop.

- Switch off the device (1.4).

- Remove the plastic cover located at the position (1.5) and pull the cord to open the door.

- Put the cord into the original position and the plastic cover again.

Don't use manual opening with the centrifuge switched on.

◆ Install the rotor. Use the designated rotors of this manual (Table 1) to avoid bad results or damage the device. Follow the instructions (fig 5):

1. Clean the surface of the motor spindle, after apply lubricant.

2. Aim the rotor central bore at the motor spindle, lower vertically until it reaches the bottom of the cone. Let your hands go and press the rotor downward.

3. Use the hexagonal 6 mm wrench provided with the device and assembly the rotor. Be sure whether installed position of rotor change or not, if necessary, re-tighten the screw to ensure solid installation of the rotor.

◆ Place the adaptor you want use and the centrifuge tubes. Verify that the whole setup can swing freely, if not could cause accidents and damage the device.



Figura colocación adaptadores

Note: The adaptors selection must be made considering its maximum rated speed or FCR values (Table 1), never exceed these values.

Note: The centrifuge tubes can load the samples at different time but it must be symmetric loading of samples (the weight of error allowed is ≤ 1.5 g). Start the machine with asymmetric loading is not allowed!

Calculation of maximum rotor load

There is huge centrifugal force when the centrifuge Nahita Multibas is running with high-speed. The design of each rotor requires sufficient mechanical strength under maximum rated speed that is "safety factor", however, this "safety factor" requires that the

load of the rotor shouldn't exceed its maximum rated load.

If the sum of sample and container exceed the maximum rated load of the rotor, you shall reduce the weight of the sample or calculate the allowed-speed of the rotor (NPERM) to make sure that the load of the rotor won't exceed its maximum rated load.

$$NPERM = N_{max} \times (\text{maximum permissible load/actual load}) \times 0.5$$

The density of the separated sample can not exceed 1.2 g/cm³, when the density of the sample exceeds 1.2 g/cm³, the maximum speed of the rotor must be reduced accordingly.

If the rotor runs unbalanced in the process operation, the spindle will overran the allowed range, and the machine will stop to run and make an alarm and the centrifuge will stop.

◆ Press the door cover pressing down until you can hear sound like "clack" (If the door is not close properly, an alarm will sound and the centrifuge will stop).

Attention! Do not start the machine with the cover door open.

◆ Parameter setting. Press the setting knob key (2.4) for set each parameter. When one parameter is selected it will blink, set the value turning the knob, after 2 or 3 seconds it will be set and the panel will come back to the start mode. The parameters that can be set are:

1. **Rotor code.** Select the rotor code for the adaptor in use (2.14) (see table 1).

IMPORTANT! ROTOR PROGRAMMED CODES MUST EXACTLY MATCH WITH THE CAPACITY OF THE TUBE AND THE ADAPTERS USED AS INDICATED IN TABLE 1

2. **Speed (RPM) or Relative centrifuge force (RCF).** Choose in section (2.9) the speed or FCR for work, selected one of the parameters is set automatically the second.

When the device is running at maximum speed may suffer a slight vibration; It is not recommend work for long periods of time at this speed, that could affect the average life of the centrifuge.

3. **Time (Hours/minutes/seconds):** There is no continuous mode operation so is necessary to program a period from 1s to 99h/59/min/59s (2.8).

4. **Setting of speed up and speed down:** Select the value between 0-9 (2.5), the bigger the value is, the shorter the time needed to get the speed.

5. **Operation mode:** Select the program (1 to 9) (2.6), the device will memorize it.

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